
DCK med



Your Specialist for Medical Diagnostic Systems



FUTREX®

6100A/ZL • 6100/XL

NIR-Technology







DCK med



DIAGNOSIS of the Body Composition with patented Near-Infrared Light-Technology

"The evidence is overwhelming that obesity affects your health and life considerably"

(NIH: National Institute of Health).

The regular DIAGNOSIS of body composition, muscles, total body water content and the basal metabolic rate (Kcal) is the ideal, and proven recommended method which has been used for years by doctors to evaluate results for oriented professional nutritional and exercise programs.

The patented FUTREX® Near-Infrared Technology has been used for over 25 years, in the field of medical technology which has the highest accuracy and repeatability results. Through fast, exact and mobile measurements, the FUTREX® diagnostic systems are especially suitable for daily professional use with your clients and patients .

DATASHEET for FUTREX® Models

| FUTREX ® Models | | FUTREX® 6100A/ZL | FUTREX® 6100/XL | FUTREX® 6050A/Pc | FUTREX® 6050/XPc |
|---|-------------------------|--|---|--|--|
| • Tech | nnology | NIR-Near Infrared Light | NIR-Near Infrared Light | NIR-Near Infrared Light | NIR-Near Infrared Light |
| • Appl | lication | Children from 5 yrs. and Adults 18 yrs. and older | Adults 18 yrs. and older | Children from 5 yrs. and Adults 18 yrs. and older | Adults 18 yrs. and older |
| | measurement uirments | No requirements for the NIR-Technology | No requirements for the NIR-Technology | No requirements for the NIR-Technology | No requirements for the NIR-Technology |
| • Mea | surement-time | appx. 2-3 min | appx. 2-3 min | appx. 2-3 min | appx. 2-3 min |
| • Cons | sumables | Thermal Paper | Thermal Paper | None | None |
| • Pow | er source | A/C Adapter/AA Battery | A/C Adapter/AA Battery | A/C Adapter | A/C Adapter |
| • Colle | eted Data | Fat Mass (kg & %) Fat Free Mass (kg & %) Total Body Water (I & %) Body Mass Index Basal Metabolic Rate (kcal) | Fat Mass (kg & %) Fat Free Mass (kg & %) Total Body Water (I & %) Body Mass Index Basal Metabolic Rate (kcal) | Fat Mass (kg & %) Fat Free Mass (kg & %) Total Body Water (I & %) Body Mass Index Basal Metabolic Rate | Fat Mass (kg & %) Fat Free Mass (kg & %) Total Body Water (I & %) Body Mass Index Basal Metabolic Rate |
| • PC - | Connection | Optional | Optional | Requirement | Requirement |
| • Soft | Software | Inclusive | Inclusive | Inclusive | Inclusive |
| Part of the state | | Fully automatic measurement wizard, assessment of body composition (total body and local analysis), plan for desired weight, illustration of the current results, plan & target, diet and exercise tips, individual printouts, a schedule for retest and follow up, Multilingual (English, Spanish, French, Italian, Dutch, Turkish, Polish, Portuguese). Additional languages can be added as requested in no time at all | | | |
| • Addi | itional Options | Data storage Programing of Company name & Advertisement Text | Data storage Programing of Company name & Advertisement Text | None | None |
| • Servi | ice | On site Training or via remote Internet access | On site Training or via remote Internet access | On site Training or via remote Internet access | On site Training or via remote Internet access |
| • Price | es | Details upon request | Details upon request | Details upon request | Details upon request |

DCK med GmbH

Echterdinger Str. 111 D-70794 Filderstadt Tel.: +49 (0)711 – 9077457-13

Tel.: +49 (0)711 – 9077457-13 Mobil: +49 (0)172 – 7291171







